



Carefully center the comfort **shield molding** with the brim of your hat, and firmly push into place.



Continue pushing the shield molding into place, following the shape of the brim and ensuring that the molding remains centered.



Repeat the same step on the other side, and ensure both sides have a snug fit with no space between the hat and comfort Shield.



Carefully peel the protective layer off both sides of shield. You are then ready for use.

## Product Advisory

### Comfort Shield's are suitable for Medical use

For multi-use wash with soap and water wash as required by the CDC and your medical facility.

Wearing our **Comfort Shield** alone will not prevent infection. It is designed to protect your eyes and face, but to prevent infection we recommend use with other personal protective equipment such as an N-95 mask.

FDA compliance under E.U.A. COVID-19 act. **Meets or exceeds ANSI Z87 BGZ87** requirements as Non-impact face shield

**GlowbackLED**®

16409 NW 8<sup>th</sup> Ave

Miami Gardens, FL 33169

1-800-679-4902 • info@bluegateinc.com

GlowbackLED® is a registered trademark of Bluegate Inc  
©2020 Bluegate Inc. All Rights Reserved  
Specifications subject to change without notice





Carefully center the comfort **shield molding** with the brim of your Hard Hat, and firmly push into place.



Continue pushing the shield molding into place on both sides, following the shape of the brim and ensuring that the molding remains centered.



After removing the protective film, your hard hat with comfort shield is ready for use.

### Using Comfort Shield With Mask Gripping Cap



**First** put your mask on. **Then**, put your **Comfort Shield with Mask Gripping Cap** on.



While wearing **Comfort Shield with Mask Gripping Cap** and face mask, pull face mask ear loops back, and attach them to the **Comfort Shield Mask Buttons**.



Ensure both mask loops are securely attached to **Comfort Shield Mask Button** on both sides of **Mask Gripping Cap**, you are then ready for use.